

Kampa watches in the background. Your tremor, sleep, and walking are tracked automatically - no tapping, no daily check-ins. The more context it has - medications, meals, workouts - the sharper its insights get, and almost all of that already lives in Apple Health and your Watch, so there's little new to keep up.

Set up — 3 steps, once

1. **Install Kampa** from the TestFlight email invitation, then open it.
2. **Allow all three permission screens** - **Health** (tap *Turn On All*), **Medications** (tap *Allow*), and **Motion & Fitness** (tap *Allow*). Motion & Fitness powers tremor tracking - without it your tremor graph stays empty while everything else works. Fix a wrong tap in *Settings > Privacy & Security > Motion & Fitness* (or *Settings > Apps > Health*).
3. **Open Kampa on your Apple Watch once** to start tremor tracking (install it from the Watch app on your iPhone if needed), then wear the watch through the day - and overnight for sleep.

Tracked automatically — you never log these

Tremor & dyskinesia from the Apple Watch; **walking & gait, sleep, heart rate, and HRV** from your iPhone and Watch - all on their own once permissions are on. Tremor starts after about a day of Watch wear, so an empty view at first is normal; if it's still empty after a day or two, check that **Motion & Fitness** is on and that you've opened Kampa on the Watch.

Wear a glucose monitor? Any brand that writes to Apple Health — Dexcom Stelo, Abbott Lingo, FreeStyle Libre — is read automatically and shown as its own band on your Review timeline. Nothing to set up beyond the CGM's own app.

The context that sharpens it — mostly things you already do

Kampa doesn't make you keep a separate logbook - it reads what you already track, so you log the way you normally do:

- **Medications** — record your doses in **Apple Health > Search > Medications**. Kampa reads them automatically - dose timing is one of its most useful insights.
- **Workouts** — just start one on your **Apple Watch** as usual; Kampa picks it up automatically.
- **Food, gut symptoms & mindfulness** — the things Health can't capture for you. Tap + on the Review screen and *say it* in plain words: "5 almonds with tea," "bloated after lunch," "ten minutes of breathing" — Kampa sorts each into the right place. Prefer typing? You can. No counting, no forms.

Reading it

- **Review** (home) — your day on one timeline: sleep, tremor, and events (💊 medication, 🍴 food, 🏃 workout). Use the ◀ ▶ arrows to change days.
- **What to expect today** — on today's Review, Kampa projects your likely *on* and *wearing-off* windows from the doses you've logged, so you can see the shape of the day ahead.
- **Insights** (💡 top-right) — patterns across weeks and months, each with an honest confidence level (*Watching* → *Emerging* → *Strong*) and a "discuss with your neurologist" note where it matters. The longer you use Kampa, the more it can tell you.

Telling us what you think

You're one of the first people to use Kampa, so your reactions matter. Send feedback from **TestFlight** (Kampa > **Send Feedback**, screenshots welcome) or just message me. No note is too small.

Requirements: iPhone 11 or later (iOS 26) and Apple Watch Series 4 or later (watchOS 10) - the Watch is how Kampa reads tremor. **Not medical advice.** Kampa shows you your own data (n-of-1); it is not a medical device and does not diagnose or treat. Never change a medication or dose based on the app - always go through your neurologist. A calm day or "no clear pattern yet" is real information, not a failure.